

# Transient Absorption Spectroscopy (TAS)

Igor Dmytruk

Department of Experimental Physics, Faculty of Physics,  
Taras Shevchenko National University of Kyiv

Department of Photonic Processes, Institute of Physics,  
National Academy of Sciences of Ukraine




# Outline

- Methods of time-resolved spectroscopy:
  - Time-resolved luminescence
  - Pump-probe:
    - Absorption (TAS)
    - Reflection
    - Anisotropy
- TAS measurement setup
  - Physical principles
  - What values are actually measured
  - Practical realization
  - Pitfalls and workarounds
- Examples of results
  - Dye molecules
  - Nanoparticles
- More sophisticated schematics


# TIME-RESOLVED SPECTROSCOPY

## Classification of Experimental Methods



# TIME-RESOLVED SPECTROSCOPY


## *Classification of Experimental Methods*




## Main idea:

EXCITATION: 150 fs, 400 nm, 1 kHz, up to  $10^9 \text{ W/cm}^2$

PROBE: 150 fs “white light”, 200 - 1600 nm



Range of delay between pump and probe: 0 fs – 2 ns



1 – pump  
2 – probe  
3 – reference

# Physical principles

- Full optical “gating”
- substitution "time -> distance“
- Absorption, transmission or reflection spectroscopy
- nonlinear element – object of study
- contrast is provided by frequency and spatial separation of pump and probe

## Data processing

measured quantity – change of optical density of the sample after pump pulse


$$\Delta D = D - D_0 =$$
$$= \log \frac{I_0 - I_{background}}{I_{pump-probe} - I_{pump}} - \log \frac{I_0 - I_{background}}{I_{probe} - I_{background}}$$

## Requirements to samples


- Stability in pump beam
- low light scattering
- relaxation time < 1ms

# Practical implementation

- Optical schematic





# Experimental setup






# Pitfalls and workarounds

- Sample instability
- Chirp of white light
- Unwanted optical signals (luminescence, scattering)
- Misalignment as delay line moves
- Nonlinear phenomena in the sample






# Femtosecond Laser Complex for collective use at Institute of Physics of National Academy of Sciences of Ukraine




# Characteristics




- time range – up to 2 ns
- repetition rate up to 1 kHz
- scanning step – 0.1  $\mu\text{m}$  ( $\sim 0.7$  fs)
- time resolution –  $\sim 150$  fs
- spectral ranges:
  - pumping – 750-800-900 nm (1 mJ)
  - 375-400-450 nm (0.2 mJ)
  - 250-267-300 nm
  - OPA and its harmonics, DFG
  - probing – 250-1100-1600 nm




25 /



# Time-resolved luminescence




Optical Kerr gate



# Electronic Transitions in a Dye Molecule Studied by TAS

*Key pathways contributing to pump–probe (transient absorption) signals*




TAS monitors changes in absorbance ( $\Delta A$ ) as a function of probe wavelength and pump–probe delay time, reflecting population dynamics among these electronic states.

# Transient Absorption Spectra and Time Traces of Dye Molecules

## Examples of TAS Measurements

### Broadband TAS Spectra ( $\Delta A$ )



### How to read TAS spectra

#### Positive $\Delta A$ (above zero):


- Excited State Absorption (ESA)
- Ground State Bleach (GSB)

#### Negative $\Delta A$ (below zero):

- Stimulated Emission (SE)



### Kinetic Time Traces ( $\Delta A$ vs pump-probe delay time)



### What time traces tell us

- Rise and decay of excited state absorption (ESA)
- Recovery of ground state (GSB)
- Stimulated emission (SE) dynamics
- Extract lifetimes and kinetic models



TAS provides a “molecular movie” of excited state processes:

✓ Excitation


✓ Relaxation

✓ Energy transfer

✓ Charge transfer




✓ Nonradiative decay

# Gain in Rhodamine 6G




# Samples preparation and characterization



- Sol-gel technique based on hydrolysis of tetraethoxysilane
- doping with Cu salt
- addition of silica powder
- gelation, drying
- annealing in air and/or hydrogen




# Absorption spectra




# Copper nanoparticles in SiO<sub>2</sub> matrix



# Copper nanoparticles in SiO2 matrix




# Colloidal solution of CdSe nanoparticles




←  $t = -3$  ps



$t = 3$  ps →




# Colloidal solution of CdSe nanoparticles




# Ultrafast spectroscopy of copper nanoparticles in silica matrix





## Kinetics of induced transparency






### Plasmon peak parameters




# Observed phenomena

- surface plasmon excitation and scattering
- band-to-band transitions ( $d \rightarrow p$ ,  $d \rightarrow s$ )
- heating and thermalization of electron gas
- electron-lattice relaxation
- excitation of vibrational mode
- energy transfer to matrix

## Time evolution of electron gas temperature




## Polarization measurement results



Plasmon peak position

Half-width of the plasmon peak




# Modern Time-Resolved Absorption Spectroscopy: Advanced Techniques

*Beyond conventional broadband TAS*

## 1. 2D-TAS (Two-Dimensional Electronic Spectroscopy)


Correlates excitation and detection wavelengths



- Resolves overlapping states and couplings
- Maps energy transfer pathways
- Reveals coherent dynamics

## 2. Femtosecond Stimulated Raman Spectroscopy (FSRS)


Vibrationally resolved TAS with high selectivity



- Mode-specific dynamics
- Probes structure and environment
- Combines electronic & vibrational information

## 3. Broadband Mid-IR TAS


Probes vibrational & charge-transfer states in the IR



- Access to vibrational fingerprints
- Sensitive to environment & H-bonding
- Tracks charge separation in the IR

## 4. High-Repetition-Rate / GHz TAS


High SNR for weak signals and long time dynamics



- Enables  $\mu$ s-ms dynamics
- Excellent SNR
- Ideal for weak absorbers and photovoltaics

## 5. Polarization-Resolved TAS (Anisotropy)


Measures rotational and energy transfer dynamics



- Rotational dynamics
- Energy transfer by dipole-dipole coupling
- Distinguishes static vs dynamic disorder

## 6. Upconversion TAS (Ultrafast & NIR)


Extends into NIR with femtosecond temporal resolution



- Sub-100 fs time resolution
- Access to NIR states
- Ideal for charge-transfer & excitons



Why these advanced TAS techniques?



Reveal couplings and energy pathways



Access vibrational, rotational & electronic degrees of freedom








Probe weak transitions with high sensitivity



Connect structure to function on ultrafast time scales

# Quinoline<sup>[1]</sup>



Thank you for attention!